

3+1



Luxator®

The Luxator product line offers a safe and precise solution for successful extractions. The fine tapered blade compresses the alveolar, cuts the membrane, and gently eases the tooth from its socket reducing damage to surrounding tissue, keeping a better anatomy for an implant site. The shape and size of the handle is designed to minimize the force needed during the extraction.



3+1

Prophylaxpasta®

Prophylaxpasta offers and delivers innovative solutions for such demands. It is one of the most well known brands worldwide, offers reliable results ranging from regular cleaning of the teeth to more complicated implant maintenance procedures.

DIRECTA

Troll® Byte Kimera

TrollByte Kimera has everything you need.

TrollByte Kimera is an aiming system for intraoral x-ray, allowing you to take any type of x-ray using only 3 holders. The one-piece design requires no assembly which facilitates its handling and it is compatible with most common sensor brands and image plates on the market. TrollByte Kimera is bio-based on renewable raw materials, and all parts are autoclavable. Autoclavable.

Colour coded for easy identification:

- Yellow - Posterior PA & Horizontal Bitewing
- Blue - Anterior PA & Vertical Bitewing
- Red - Horizontal Bitewing



Arm

The stable arm makes it easy to keep the sensor in place.

Aiming Ring

The aiming ring is colour coded to match the holders and has 5 different image positions for easy use. You only need one ring for all holders.

Claw

The size of the claw is designed to fit most sensor brands on the market, with the possibility for small adjustments when needed. The claw has round edges for patient comfort.

Troll® Foil

TrollFoil takes the guesswork out of occlusal adjustments.

TrollFoil can be used under a wide variety of clinical situations including wet or dry teeth, limited opening, limited vestibular space, gaggers, and metal and non-metallic restorations. You are able to verify occlusal contacts for both arches from one tooth to an entire quadrant using the mini or the original quadrant size. The double-sided foil is only 8 microns thick, and it has no problem marking wet surfaces, dry surfaces or highly polished surfaces such as cast gold or BruxZir.




TrollFoil campaign

The only articulating foil you will ever need!



DIRECTA



20%