




Introduction

Simplicity & Versatility

Practice Growth

Patient Care

Peace of Mind




Primescan[®] 2

Enter a new era of digital patient care

Introducing Primescan 2
Powered by DS Core[®]

Our direct-to-cloud, wireless intraoral scanner that stands for simplicity and versatility, helps you grow your practice and deliver excellent patient care.





Take your next step in digital dentistry with Primescan 2

The first cloud-native intraoral scanning solution



Primescan 2 is an easy way to evolve your practice, no matter where you are in your digital journey.



Your versatile day-to-day companion

- Easy to integrate into your existing set-up to scan on any computer in your practice
- With a simple direct-to-cloud scanning workflow that quickly releases Primescan 2 after each scan



Helps grow your practice

- With access to diverse indications, so you can offer more to your patients
- Enabling effective patient communication via DS Core¹ to support you with treatment acceptance



Supports you in delivering exceptional patient care


- With the speed and accuracy of the proven Primescan technology, now wireless and redesigned for more patient comfort
- Caries detection aid to support you with timely diagnosis²

¹ Primescan 2 requires a paid subscription to DS Core.

² Currently not available in the CE markets, UK and Switzerland. Please refer to your Dentsply Sirona representative to confirm availability in your region.



Scan any patient, anywhere, anytime



Beyond wireless

Primescan 2 is powered by DS Core. The scan data is captured directly on the cloud platform, giving you the flexibility to use your scanner anywhere, anytime you need it.

Scan anywhere you need

Primescan 2 doesn't require any dedicated hardware - use it with any computer, laptop, tablet or other mobile device with an internet connection. Effortlessly transform every room into a scanning room.

One scanner for numerous treatments

Primescan 2 delivers excellent accuracy and reliable scan quality no matter the indication or material. Thanks to its broad range of applications, it can open up a world of possibilities for your practice offering.

Instantly ready to scan again

As soon as the scan is complete, Primescan 2 can be moved to the next patient, allowing for a fast hand-over of the scanner within the practice and increasing your device's availability.

"Primescan 2 is way more than a wireless scanner, it's all I need to scan every patient. Not only can I use it anywhere in my practice without being tied to a specific computer, but I can also confidently apply it to a wide range of treatments, from basic restorative cases to more complex interventions and specialized procedures."




Dr. Marko Ahonen³ | Dentist | Finland

³The professionals giving testimonials on this document may have been compensated for use of their experiences.



"Primescan 2 allows me to start every patient consultation with a scan. I can take a status scan at first or a treatment scan if treatment has already been determined. This helps me **save time**, quickly frees-up my scanner, allows for a shorter scan experience for the patient, and enables me to **easily use the scan data across cases.**"


Taking intraoral scanning to a whole new level with **DS Core⁴**



Start every patient consultation with a scan

Take a status scan to check-up your patient's situation, then take a treatment scan to focus on a specific area to treat.


This flexibility helps you to reduce chair time and lets you use scan data across cases in DS Core.



Easy collaboration with lab partners

Once you finish the scan, easily order with just a few clicks, share files, send messages and signed lab proposals in real-time for more effective communication. Or have the practice team prepare the order draft for your review.


The advantage of using any device to collaborate with your lab gives you the chance to better delegate to your team.



Efficient patient communication

Visualization with DS Core helps your patients understand their treatment need.


On Canvas, you can see, arrange and annotate scans along with any other images in a single comprehensive view. Plus, you can easily share the files with patients for their review anywhere.



One scan – full flexibility

Launch your next workflow flexibly from one space – Simulate an aligner treatment, send to your lab, open in CEREC or transfer to your specialist.

Primescan 2 treatment scans automatically create a treatment space to accelerate your next step and allow you to manage all treatment related files in one central space.



All without the headache of software updates

As a cloud-based platform, DS Core ensures you're always using the latest software version every time you log in, without needing to install updates or upgrades.

We are continuously rolling out innovative features to streamline and improve your workflow experience.



Dr. Verena Freier⁵
Dentist | Germany

⁴ DS Core is not a medical device software. It is not intended to diagnose, treat, cure, or prevent any disease or health condition. Use only for information, education or sharing purposes. Service is subject to country availability. Please contact your local DS representative for more details.


⁵ The professionals giving testimonials on this document may have been compensated for use of their experiences.



Unlock new treatment opportunities with Primescan 2 and take your practice to the next level

Primescan 2 together with DS Core is your easy starting point to help to grow your practice.

You can add new treatments to your portfolio at your pace, including the opportunity to expand into in-house manufacturing.




COLLABORATION

Labs & Partners

Work with any lab or partner



- Indirect restorations
- Implant restorations



- Clear aligners



- Implantology
- Guided surgery



- Sleep devices
- Splints



- Dentures



Seamless order with SureSmile® or work with any lab



Atlantis® Siplant®



Panthera™ ResMed




Reference dentures


CHAIRSIDE

In-house Manufacturing

Fabricate your own restorations and appliances



- Indirect restorations
- Implant restorations



- Surgical guides, models, night guards, temporary restorations & dentures



CEREC® inLab®






Primeprint®

Click on the treatment icons to learn more about each workflow



Restorative workflow with Primescan 2

Primescan 2 is an excellent choice to bring more digitalization, predictability and efficiency in your restorative workflow. You and your favorite dental lab partner benefit from highly accurate scans, easy data transfer and simple collaboration on DS Core. Continue what you're doing today - just better.



Scan

Scan your patient with Primescan 2 and the intraoral scan is captured in DS Core - freeing up your scanner immediately.

Order

Place an order with your lab of choice from wherever you are via DS Core.

Easily collaborate with your lab

Keep track of your orders, give guidance to your lab via chat, approve designs or add files if needed. Your lab can easily open files with DS Core and work with the design software of their choice.

Treating your patient

With Primescan 2, you can provide a broad variety of restorative indications giving your patients the treatment they deserve, offering them indirect restorations and esthetic treatment options.



81%

of Primescan clinicians agree that Primescan improves the quality of final restorations.⁶



Clear aligner workflow with Primescan 2

Leverage Primescan 2 for a streamlined process with SureSmile® Clear Aligners, enabling faster aligner ordering. Enhance patient engagement and help drive treatment plan acceptance with the SureSmile® Simulator's 3D visualization.



Scan

Scan your patient with Primescan 2 and the intraoral scan is captured in DS Core - freeing up your scanner immediately.

Simulate & Show

Use the SureSmile® Simulator⁷ to help increase chairside acceptance of aligner treatment by creating a 3D-visualization of your patient's possible new smile.

Share the simulation with your patient for them to review at home.

Submit

Start the order for your patient's aligner treatment plan including all the collected records stored in DS Core.

Review & Approve

Review the SureSmile treatment plan on DS Core to easily approve and finalize the order for physical aligners.

For modification, you can easily navigate back to SureSmile to utilize the hands on tools you love.

Start Treatment

When the SureSmile® Aligners, Patient Kit and optional SureSmile® VPro™⁸ arrive at your practice, it's time to start treatment.

81%

ofsmile-conscious people say they'd like their teeth to be straighter - indicating many patients in your dental practice could benefit from clear aligner therapy.⁹

⁷ Requires user to have a paid DS Core subscription and Primescan. DS Core only works with Sure Smile aligner software platform. Displayed results are a simulation of potential aligner prescribed treatment plan, which must be directed by the dentist and complied with by the patient. Actual treatment outcomes may differ substantially.




⁸ Not available in every country. Talk to your sales representative for market availability.

⁹Data on file. Available at contact@dentsplysirona.com



Implantology workflow with Primescan 2

Start your implant workflow with Primescan 2 and benefit from highly accurate scan data and effective communication with your lab. Directly connect to great partners like Simplant® and Atlantis®. Make use of various workflow options for both implant surgery and implant prosthetics.



Scan

Simulate & Show

Bring all files together

Collaborate

Treating your patient

Scan your patient with Primescan 2 and the intraoral scan is captured in DS Core - freeing up your scanner immediately.

Use any device with a web browser to access and arrange relevant patient files for a holistic consultation on DS Core. Help your patients understand the treatment plan by highlighting findings and communicating the treatment need with added annotations.

Bring all the relevant case files together in your patient media library. Your intraoral scans with Primescan 2 are available immediately in the patient media library and you can add other relevant case files like X-ray files.

Collaborate with labs or refer patients to an implant specialist by easily sharing patient files. The two-way media sharing and chat functionality allows you to get feedback efficiently.

Benefit from having a high-level of control over your workflow while saving time and providing your patients with treatment comfort.

"Primescan 2 gives me an **unmatched level of accuracy** that I truly appreciate with all of my restorations. The scans are remarkably clear and precise thanks to its dynamic depth scan technology, at a scan depth of up to 20 mm deep, it becomes a **must have** with implant restorations. No matter how complex the treatment, I can be **confident** that the **results will be excellent**, and my patients will receive the best care I can offer."



Dr. Daniel Butterman¹⁰ | Dentist | USA

¹⁰ The professionals giving testimonials on this document may have been compensated for use of their experiences.



Integrate Primescan 2 into your CEREC workflow



Whether you are already working with CEREC or you would like to do it in the future, Primescan 2 can expand your offering for restorative treatments and allows you to turn every room in your practice into a CEREC room. So you can scan anytime, anywhere you need it.

With Primescan 2 you have the flexibility to expand into in-house manufacturing when you are ready. The scan can be opened from DS Core in CEREC Software and integrated into your workflows of tomorrow.

By adding Primescan 2 to your CEREC workflow, you have the option to continue designing restorations on your Primescan AC, or utilize a design station with DS Core. This ensures that the scanner is immediately ready for the next patient.





Scan

Scan your patient with Primescan 2 and the intraoral scan is captured in DS Core.



Free up scanner for next patient

Initiate your next workflow step in DS Core from any device - without blocking the scanner.



Design on your preferred design station

Design the restoration using CEREC SW 5.3, or inLab SW 22.3, either on your Primescan AC or on a dedicated design station.



Manage your practice workflows

Efficiently organize and manage your practice workflows with DS Core and access files from any device with a browser. Experience the advantages of simplified equipment and manufacturing management¹¹.

"Primescan 2 has turned every treatment room in my office into a functional scanning and CEREC room. It allows me to wirelessly scan all my patients with ease directly to the cloud for patient records, treatment planning, or designing a restoration on my multiple workstations. It has truly changed the way I practice resulting in excellent patient flow and overall efficiency in my office."



Dr. Mike Skramstad¹² | Dentist | USA

¹¹ Manufacturing on DS Core requires a paid subscription to DS Core.

¹² The professionals giving testimonials on this document may have been compensated for use of their experiences.



Designed to be fun for you and your team



Making the day-to-day fun

Primescan 2 is an easy way to bring innovation into your practice. The ergonomic and intuitive design enables easy, fun handling - every day.



Intuitive to use, quick to learn

Primescan 2 is easy to learn, so your practice team can scan and work with files in no time.



Empower your practice team

With Primescan 2 you can easily delegate scanning and ordering to your staff, so you can focus on what matters the most - your patients.

Your practice team will appreciate working with Primescan 2 and feel confident with the use of a state-of-the-art technology, enriching their work routines. Empowered team members stay ahead.



Keep your team engaged

According to a recent study¹³, dental practice staff find digital equipment reassures that they are giving the best patient care possible and prefer to work in a digital practice.

¹³Data on file. Available at contact@dentsplysirona.com

¹⁴ The professionals giving testimonials on this document may have been compensated for use of their experiences.

"As a hygienist, taking digital impressions with Primescan 2 for our patients wellness checks makes my **daily work very dynamic**. It's very **intuitive** and comfortable to handle. I can change catalogues by just tapping anywhere on the scanner. The greatest benefit is that it aids with the **work-flow between the doctor and I** because they can use the scans for treatment as well. I feel proud and happy to be able to provide excellent care to my patients by using innovative technology."



Stephanie Lucia Vasquez ¹⁴ |
Dental Assistant | USA



Enhance your patients' experience

80%

of patients feel more confident that they are getting excellent dental care when you choose digital technology such as an intraoral scanner.¹⁵

84%

of patients prefer digital to conventional impression taking.¹⁶




Minimize scanning time

Primescan 2 completes a full arch scan in less than 1 minute, so your patients don't have to spend a long time in the chair being scanned.¹⁷



Great patient comfort

The size and outer shape of our scanner tip has been redesigned with a leaner angle profile to easily scan the molars and hard to reach distal surfaces, making a comfortable scanning experience for your patients.



Caries detection and prevention with timely diagnose

Early diagnosis is essential for maintaining oral health. Using fluorescence and near-infrared, Primescan 2 can assist with caries identification¹⁸, while visualization with DS Core helps your patients to understand their treatment need¹⁹.

¹⁵ US Patient Survey, May 2020, N=254 US Adult Dental Patients.


¹⁶ Haddadi Y, Bahrami G, Isidor F, Evaluation of Operating Time and Patient Perception Using Conventional Impression Taking and Intraoral Scanning for Crown Manufacture: A Split-mouth, Randomized Clinical Study. Int J Prosthodont 2018;31(31):55-59

^{17,18} Data on file. Available at contact@dentsplysirona.com

¹⁹ Currently not available in the CE markets, UK and Switzerland.

Please refer to your Dentsply Sirona representative to confirm availability in your region.

Caries Detection Tool






Proven scanning technology for high quality patient care

- Expect excellent clinical results thanks to high-quality, accurate scans from the renowned Primescan technology²⁰
- The new ergonomic design of Primescan 2 makes it easy to scan difficult-to-access areas or patients who have a small oral cavity²¹
- Provide excellent patient comfort with Primescan 2: you can complete a full arch scan in less than a minute²²

Click on the light blue icons to find out more about the established Primescan technology.

Click on the dark blue icons to find out more about the new technological features of Primescan 2.



²⁰ Primescans' precision and accuracy have been proven in several studies. For more detail, please refer to the [Dentsply Sirona Clinical Study Overview](#)

^{21, 22, 23, 25} Data on file. Available at contact@dentsplysirona.com

²⁴ Currently not available in the CE markets, UK and Switzerland. Please refer to your Dentsply Sirona representative to confirm availability in your region.



Comprehensive hygiene concept

The closed design around the field of view and sleeve concept leave no doubt about fulfilling hygiene needs and complying to regulatory requirements. Protect your patients from contamination risk with a comprehensive hygiene concept with different sleeves that fit your practice needs.

Designed for your patient's comfort

The size and outer shape of the scanner tip have been redesigned to make scanning comfortable, especially for patients with a smaller mouth or those having difficulty to open the jaw. Plus, the smaller tip size makes it easier to reach molars and distal surfaces.

²⁵ Primescan 2 Multiuse Sleeve available in second half 2024.


²⁶ Please consider country specific regulations for hygiene.



The **single-use sleeve** helps to make intraoral scanning very comfortable and safe. Simply change the sleeve after each patient scan.



The **Primescan 2 Multiuse Sleeve^{25,26}** is a steel sleeve with sapphire glass window that can be treated with High Level Disinfection (HLD).





Feel reassured with DS Core Care



DS Core Care is the advanced technical service and support solution to help your equipment and practice run reliably and efficiently.

Always in control of your equipment.

By protecting your equipment, DS Core Care helps you to optimize your Primescan 2 availability and uplift your productivity.

- Express replacement
- Spare parts
- Accidental damage coverage
- Hotline support

Thanks to the connection to DS Core, you receive guided support on DS Core to get your Primescan 2 ready for operation - and real-time troubleshooting for your set up questions. So, you and your team feel confident using your new intraoral scanner.

- Fast Service Clarification
- Customer Support Portal



Individualize your Primescan 2 experience

Primescan 2 accessories allow you to configure your scanner to your specific needs.

Primescan 2 Cart



 **Ergonomic design**
A slim and versatile design for smooth mobility around your practice.


 **Easy patient communication**
An all-in one solution always in direct reach for effortless and effective consultations.

 **Always scan-ready**
Wireless connection and charger for a simple scan experience.

Exchangeable battery concept

With rechargeable batteries, use Primescan 2 for up to four scans per charge (per battery)²⁷

- Charging state display
- Batteries are also available for wall mount and Primescan 2 Cart
- One charger and 3 batteries included - additional chargers can be ordered separately.



Dedicated cradle for stability and care

- Primescan 2 comes with a cradle that safely holds the scanner while not in use
- For additional care of your scanner, you can place a cradle in every room where you will use Primescan 2

²⁷ Data on file. Available at contact@dentsplysirona.com



Are you ready to enter a new era of digital patient care?

For more product information or to schedule a Primescan 2 demonstration, please contact your local Dentsply Sirona representative.



Dentsply Sirona


SIRONA Dental Systems GmbH
Fabrikstraße 31, 64625 Bensheim,
Deutschland
dentsplysirona.com

Subject to technical changes and errors in the text. Order No. A91100-M43-C372-01-7600.
Registered trademarks, trade names and logos are used. Even if these are not identified as such in the respective places, the corresponding legal provisions apply. Unless otherwise stated, all comparative statements in this document refer to a comparison of Dentsply Sirona products with each other.



Proven scanning technology for high quality patient care

- Expect excellent clinical results thanks to high-quality, accurate scans from the renowned Primescan technology²⁰
- The new ergonomic design of Primescan 2 makes it easy to scan difficult-to-access areas or patients who have a small oral cavity²¹
- Provide excellent patient comfort with Primescan 2: you can complete a full arch scan in less than a minute²²



²⁰ Primescans' precision and accuracy have been proven in several studies. For more detail, please refer to the [Dentsply Sirona Clinical Study Overview](#).


^{21,22,23,25} Data on file. Available at contact@dentsplysirona.com

²⁴ Currently not available in the CE markets, UK and Switzerland. Please refer to your Dentsply Sirona representative to confirm availability in your region.



Proven scanning technology for high quality patient care

- Expect excellent clinical results thanks to high-quality, accurate scans from the renowned Primescan technology²⁰
- The new ergonomic design of Primescan 2 makes it easy to scan difficult-to-access areas or patients who have a small oral cavity²¹
- Provide excellent patient comfort with Primescan 2: you can complete a full arch scan in less than a minute²¹



²⁰ Primescans' precision and accuracy have been proven in several studies. For more detail, please refer to the [Dentsply Sirona Clinical Study Overview](#).


^{21,22,23,25} Data on file. Available at contact@dentsplysirona.com

²⁴ Currently not available in the CE markets, UK and Switzerland. Please refer to your Dentsply Sirona representative to confirm availability in your region.



Proven scanning technology for high quality patient care

- Expect excellent clinical results thanks to high-quality, accurate scans from the renowned Primescan technology²⁰
- The new ergonomic design of Primescan 2 makes it easy to scan difficult-to-access areas or patients who have a small oral cavity²¹
- Provide excellent patient comfort with Primescan 2: you can complete a full arch scan in less than a minute²²



²⁰ Primescans' precision and accuracy have been proven in several studies. For more detail, please refer to the [Dentsply Sirona Clinical Study Overview](#).


^{21,22,23,25} Data on file. Available at contact@dentsplysirona.com

²⁴ Currently not available in the CE markets, UK and Switzerland. Please refer to your Dentsply Sirona representative to confirm availability in your region.



Proven scanning technology for high quality patient care

- Expect excellent clinical results thanks to high-quality, accurate scans from the renowned Primescan technology²⁰
- The new ergonomic design of Primescan 2 makes it easy to scan difficult-to-access areas or patients who have a small oral cavity²¹
- Provide excellent patient comfort with Primescan 2: you can complete a full arch scan in less than a minute²²



✕ Battery

Primescan 2 comes with 3 rechargeable batteries each can last up to 4 scans so you can scan all day.²⁰

²⁰ Primescans' precision and accuracy have been proven in several studies. For more detail, please refer to the [Dentsply Sirona Clinical Study Overview](#).


^{21,22,23,25} Data on file. Available at contact@dentsplysirona.com

²⁴ Currently not available in the CE markets, UK and Switzerland. Please refer to your Dentsply Sirona representative to confirm availability in your region.



Proven scanning technology for high quality patient care

- Expect excellent clinical results thanks to high-quality, accurate scans from the renowned Primescan technology²⁰
- The new ergonomic design of Primescan 2 makes it easy to scan difficult-to-access areas or patients who have a small oral cavity²¹
- Provide excellent patient comfort with Primescan 2: you can complete a full arch scan in less than a minute²²



²⁰ Primescans' precision and accuracy have been proven in several studies. For more detail, please refer to the [Dentsply Sirona Clinical Study Overview](#).


^{21,22,23,25} Data on file. Available at contact@dentsplysirona.com

²⁴ Currently not available in the CE markets, UK and Switzerland. Please refer to your Dentsply Sirona representative to confirm availability in your region.



Proven scanning technology for high quality patient care

- Expect excellent clinical results thanks to high-quality, accurate scans from the renowned Primescan technology²⁰
- The new ergonomic design of Primescan 2 makes it easy to scan difficult-to-access areas or patients who have a small oral cavity²¹
- Provide excellent patient comfort with Primescan 2: you can complete a full arch scan in less than a minute²²



²⁰ Primescans' precision and accuracy have been proven in several studies. For more detail, please refer to the [Dentsply Sirona Clinical Study Overview](#).


^{21,22,23,25} Data on file. Available at contact@dentsplysirona.com

²⁴ Currently not available in the CE markets, UK and Switzerland. Please refer to your Dentsply Sirona representative to confirm availability in your region.



Proven scanning technology for high quality patient care

- Expect excellent clinical results thanks to high-quality, accurate scans from the renowned Primescan technology²⁰
- The new ergonomic design of Primescan 2 makes it easy to scan difficult-to-access areas or patients who have a small oral cavity²¹
- Provide excellent patient comfort with Primescan 2: you can complete a full arch scan in less than a minute²²




²⁰ Primescans' precision and accuracy have been proven in several studies. For more detail, please refer to the [Dentsply Sirona Clinical Study Overview](#).
²⁴ Currently not available in the CE markets, UK and Switzerland. Please refer to your Dentsply Sirona representative to confirm availability in your region.



Proven scanning technology for high quality patient care

- Expect excellent clinical results thanks to high-quality, accurate scans from the renowned Primescan technology²⁰
- The new ergonomic design of Primescan 2 makes it easy to scan difficult-to-access areas or patients who have a small oral cavity²¹
- Provide excellent patient comfort with Primescan 2: you can complete a full arch scan in less than a minute²²



²⁰ Primescans' precision and accuracy have been proven in several studies. For more detail, please refer to the [Dentsply Sirona Clinical Study Overview](#).


^{21,22,23,25} Data on file. Available at contact@dentsplysirona.com

²⁴ Currently not available in the CE markets, UK and Switzerland. Please refer to your Dentsply Sirona representative to confirm availability in your region.



Proven scanning technology for high quality patient care

- Expect excellent clinical results thanks to high-quality, accurate scans from the renowned Primescan technology²⁰
- The new ergonomic design of Primescan 2 makes it easy to scan difficult-to-access areas or patients who have a small oral cavity²¹
- Provide excellent patient comfort with Primescan 2: you can complete a full arch scan in less than a minute²²



²⁰ Primescans' precision and accuracy have been proven in several studies. For more detail, please refer to the [Dentsply Sirona Clinical Study Overview](#).


^{21,22,23,25} Data on file. Available at contact@dentsplysirona.com

²⁴ Currently not available in the CE markets, UK and Switzerland. Please refer to your Dentsply Sirona representative to confirm availability in your region.



Proven scanning technology for high quality patient care

- Expect excellent clinical results thanks to high-quality, accurate scans from the renowned Primescan technology²⁰
- The new ergonomic design of Primescan 2 makes it easy to scan difficult-to-access areas or patients who have a small oral cavity²¹
- Provide excellent patient comfort with Primescan 2: you can complete a full arch scan in less than a minute²²



²⁰ Primescans' precision and accuracy have been proven in several studies. For more detail, please refer to the [Dentsply Sirona Clinical Study Overview](#).


^{21,22,23,25} Data on file. Available at contact@dentsplysirona.com

²⁴ Currently not available in the CE markets, UK and Switzerland. Please refer to your Dentsply Sirona representative to confirm availability in your region.



Proven scanning technology for high quality patient care

- Expect excellent clinical results thanks to high-quality, accurate scans from the renowned Primescan technology²⁰
- The new ergonomic design of Primescan 2 makes it easy to scan difficult-to-access areas or patients who have a small oral cavity²¹
- Provide excellent patient comfort with Primescan 2: you can complete a full arch scan in less than a minute²²



²⁰ Primescans' precision and accuracy have been proven in several studies. For more detail, please refer to the [Dentsply Sirona Clinical Study Overview](#).


^{21,22,23,25} Data on file. Available at contact@dentsplysirona.com

²⁴ Currently not available in the CE markets, UK and Switzerland. Please refer to your Dentsply Sirona representative to confirm availability in your region.



Proven scanning technology for high quality patient care

- Expect excellent clinical results thanks to high-quality, accurate scans from the renowned Primescan technology²⁰
- The new ergonomic design of Primescan 2 makes it easy to scan difficult-to-access areas or patients who have a small oral cavity²¹
- Provide excellent patient comfort with Primescan 2: you can complete a full arch scan in less than a minute²²



²⁰ Primescans' precision and accuracy have been proven in several studies. For more detail, please refer to the [Dentsply Sirona Clinical Study Overview](#).


^{21,22,23,25} Data on file. Available at contact@dentsplysirona.com

²⁴ Currently not available in the CE markets, UK and Switzerland. Please refer to your Dentsply Sirona representative to confirm availability in your region.



Proven scanning technology for high quality patient care

- Expect excellent clinical results thanks to high-quality, accurate scans from the renowned Primescan technology²⁰
- The new ergonomic design of Primescan 2 makes it easy to scan difficult-to-access areas or patients who have a small oral cavity²¹
- Provide excellent patient comfort with Primescan 2: you can complete a full arch scan in less than a minute²²



✕ Motion Sense

The scanner is switched on automatically when removed from its cradle. Scanning can start less than 1 second after automatic activation.

²⁰ Primescans' precision and accuracy have been proven in several studies. For more detail, please refer to the [Dentsply Sirona Clinical Study Overview](#).
^{21,22,23,25} Data on file. Available at contact@dentsplysirona.com
²⁴ Currently not available in the CE markets, UK and Switzerland. Please refer to your Dentsply Sirona representative to confirm availability in your region.