

Next-Generation Single-Step Adhesive

CLEARFIL™ Universal Bond Quick 2



CLEARFIL™ Universal Bond Quick 2 is the latest evolution in dental bonding technology, offering a streamlined, efficient solution for dental practitioners. Designed to enhance your workflow, this adhesive provides robust performance across a variety of clinical situations.

Choose CLEARFIL™ Universal Bond Quick 2

- ✓ Easy application
- ✓ Quick and straightforward procedure
- ✓ consistently durable bond
- ✓ high strength bonding layer
- ✓ Versatility for various applications
- ✓ Minimized risk of pooling
- ✓ Ready to use anytime

The complex technology behind simplicity

At Kuraray Noritake Dental Inc., we have refined our monomer technology to deliver exceptional bonding performance with minimal effort. Our proprietary Advanced Rapid Bond Technology combines three key monomers – the original MDP, Amide and Urethane Tetra Methacrylate. This creates a powerful adhesive that supports both ease of use and long-lasting results.



Article nr	Product	Contents
#4042-EU	CLEARFIL™ Universal Bond Quick 2 Bottle Refill	5.0 ml x 1 bottle
#4044-EU	CLEARFIL™ Universal Bond Quick 2 Bottle Value Pack	5.0 ml x 3 bottles
#4045-EU	CLEARFIL™ Universal Bond Quick 2 Bottle Trial	1.0 ml x 1 bottle
#4046-EU	CLEARFIL™ Universal Bond Quick 2 Unit Dose Standard Pack	1 ml tip x 50 tips, Applicator brush fine silver x 50pcs
#4048-EU	CLEARFIL™ Universal Bond Quick 2 Unit Dose Trial	0.1 ml tip x 5 tips, Applicator brush fine silver x 5pcs

