

Universal. Flowable. Strong. Your new go-to composite.

CLEARFIL MAJESTY™ ES Flow Universal



CLEARFIL MAJESTY™ ES Flow Universal represents a fundamental shift in the idea of what universal composites can do. By integrating high filler load and achieving high flexural strength, CLEARFIL MAJESTY™ ES Flow Universal positions itself capable of delivering aesthetic excellence as a permanent solution. Its unique light diffusion technology allows working with only one shade-U for almost all posterior and anterior cavities and shades. UD (Universal dark) shows best results in cavities of shade B3 and beyond that are lacking a back wall, such as Class IV cavities. For preferred handling, the product offers two flowability options, Low and Super low.

1 SHADE

Covers all 16 VITA* shades in all cavities except those lacking of a back wall and of shade B3 and darker.

2 FLOWABILITIES

Low and super low for a wide range of indications.


IMPRESSIVE STRENGTH

Sufficient for use even in stress-bearing posterior restorations.



All restorations. Only two shades. How does it work?

A smart combination of an enamel-like translucency, dentin-like chroma and Light-Diffusion Technology (LDT) does the trick. Thanks to the enamel-like translucency, light is transmitted through the material, bringing out the colour of the surrounding tooth structure. Fine-tuning of the chroma, which is not affected by the surrounding tooth structure, supports natural outcomes. Finally, LDT makes restorations diffuse light in a similar way as enamel and dentin. As a consequence, the material integrates easily and smoothly into the surrounding tooth structure.



| Article nr | Product | Contents |
|------------|--|-------------------------------------|
| #4600-EU | CLEARFIL MAJESTY™ ES Flow Universal Low (U) Syringe | Syringe 2.7g/1.5ml, 15x Needle Tips |
| #4601-EU | CLEARFIL MAJESTY™ ES Flow Universal Low (UD) Syringe | Syringe 2.7g/1.5ml, 15x Needle Tips |
| #4620-EU | CLEARFIL MAJESTY™ ES Flow Universal Super Low (U) Syringe | Syringe 2.7g/1.5ml, 15x Needle Tips |
| #4621-EU | CLEARFIL MAJESTY™ ES Flow Universal Super Low (UD) Syringe | Syringe 2.7g/1.5ml, 15x Needle Tips |
| #3331-EU | CLEARFIL MAJESTY™ ES Flow Needle Tip (N) | 20x |

*VITA is a trademark of VITA Zahnfabrik, Bad Sackingen, Germany

