



DS Core®

Unlocking connected dentistry
through cloud-powered workflows

✓ GROWTH


Grow your dental practice with enhanced patient communication and stay flexible with solutions that adapt to your needs.

✓ EFFICIENCY

Increase efficiency with time-saving, yet flexible cloud-powered and connected workflows.

✓ OUTCOME

Get excellent clinical outcomes and deliver a great patient experience with organizational and clinical tools for all your clinical jobs.



[Learn More →](#)

 **Dentsply
Sirona**



● Diagnose¹

● Plan

● Make

Canvas

Enhance communication and engagement

Diagnostic tools

Take CBCT length measurements with ease

Smart View

CBCT-based illustrations of your patient's anatomy²

SureSmile Simulator

Illustrate the value of a SureSmile[®] aligner treatment proposal

CEREC connected with DS Core

Initiate CEREC[®] design workflow and transfer data seamlessly⁴

CEREC on DS Core

Cloud-powered design of crowns, inlays and onlays

Manufacturing

Manage all your manufacturing jobs in a single location⁵

Ordering

Manage all your lab and service orders in a single location

IO scan & x-ray acquisition

Direct-to-cloud scanning³, connected IOS and x-ray upload

Patient & Dentist Sharing

Share media with patients and collaborate with peers

Online Media Library

Store⁶ and view patient media in one place

DS Core[®] is a Software as a Service

3m trial option⁷

Light

\$19,95

per practice

- ✓ Everything in Access
- ✓ 1 TB Online Media Library
- ✓ Connected^{8,9,10} & direct-to-cloud¹¹ scan acquisition
- ✓ Connected X-ray upload & historic sync with Sidesis⁸
- ✓ Direct PMS Sync¹²
- ✓ Patient & Dentist Sharing
- ✓ Canvas
- ✓ SureSmile Simulator
- ✓ Connected CEREC designs¹³

Standard

\$99,95

per practice

- ✓ Everything in Light
- ✓ 6 TB Online Media Library
- ✓ Primary Storage for x-rays
- ✓ Diagnostic tools
- ✓ Smart View
- ✓ CEREC on DS Core designs (crown, inlay, onlay)
- ✓ Manufacturing¹⁴ (up to 2 manufacturing devices)

Advanced

\$249,95

per practice

- ✓ Everything in Standard
- ✓ 15 TB Online Media Library
- ✓ CEREC on DS Core for Trios scans & STL export
- ✓ Manufacturing¹⁴ (unlimited manufacturing devices)

CCC replacement

Access

\$0

per practice

- ✓ Temporary storage¹⁵
- ✓ Connected scan upload with Primescan¹⁰
- ✓ Connected X-ray upload with Sidesis⁸
- ✓ Ordering

Sign up on DSCore.com →

¹ Diagnose feature availability depends on local medical device clearance; ² Smart View is not a medical device, and it is intended for patient communication only; ³ Available only with Primescan 2; ⁴ CEREC 5.3.1 license required; ⁵ inLab CAM connected to DS Core required; ⁶ DS Core offers primary storage for imaging files and Primescan 2; ⁷ You can try all the features of DS Core Light for 3 months prior to your first payment being due. You will be asked to provide important feedback about your reasons for cancelling your subscription; ⁸ Automatic upload works with Connect/CEREC SW 5.2.3 or higher and Sidesis 4 V4.3.1 or higher. Use Connect 5.2.9, CEREC SW 5.3.3 and Sidesis 4 V4.4, for best performance and upload possibilities; ⁹ Seamless upload of IO scans from 3Shape Trios 3,4,5 scanners using 3Shape Unite. Requires a DS Core paid subscription; ¹⁰ Also works with Omnicam; ¹¹ Requires Primescan 2; ¹² Requires compatible Practice Management System; ¹³ Works with CEREC 5.3.1 and higher; ¹⁴ Requires inLab CAM SW Connected with DS Core; ¹⁵ DS Core Access only provides temporary file storage. Uploaded files will be deleted after 7 days, if not used within an order. All files included in an order will be deleted 30 days once order is completed.