



# Intraoral Imaging Accessories

## Product Catalog

# Overview

**Rinn® XCP-ORA®**  
One Ring & Arm Positioning System ..... 3

**Digital Sensors**

**Universal & Autoclavable Sensor Holders**  
Rinn® XCP-DS Fit® ..... 4, 5  
Xios Intraoral Holder Systems ..... 6, 7  
Rinn® Snap-A-Ray® DS ..... 8

**Universal & Disposable Sensor Holders**  
Rinn® Uni-Grip® 360 ..... 9

**Film**

Schick AimRight Holders vs.  
Rinn XCP Holders Color Guide ..... 7

**Universal & Autoclavable Film Holders**  
Rinn® Snap-A-Ray® Xtra ..... 12  
Rinn® XCP ..... 12

**Rinn User Guide** ..... 14

**Rinn Color Guide** ..... 15

## Phosphor Plates

**Universal & Autoclavable Phosphor Plate Holders**  
Rinn® XCP-PSP Fit® ..... 10  
Phosphor Plate Envelopes ..... 11

Decrease inventory clutter  
with only **1 Arm and 1 Ring.**




## Rinn® XCP-ORA® One Ring & Arm Positioning System

**One Ring. One Arm. One Goal:  
Faster & Easier Positioning.**


- 3-in-1 positioning convenience for anterior, posterior, and bitewing
- Reduces components needed to position radiographic media
- To be used with Rinn® biteblocks for
  - Digital Sensors (XCP-DS Fit®)
  - Phosphor Plates (XCP-PSP Fit®)
  - Film (XCP®)
- Works with round and rectangular collimators
- Autoclavable

ORA Ring & Arm	Qty
550771 XCP-ORA® Arm & Ring	1 pc.
550772 XCP-ORA® Arm & Ring	2 pcs.
550773 XCP-ORA® Ring	1 pc.
550774 XCP-ORA® Arm	1 pc.

**For Digital Sensors:**



**For Phosphor Plates:**




**For Film:**




Please note: The green Endo-Biteblocks of the product lines Rinn® XCP®, XCP-PSP Fit™ and XCP-DS Fit® need the XCP® positioning arm and ring (see page 5, 9 and 11).





Easy. Hygienic. Efficient.




Xios Anterior holder



Xios Posterior holder



Xios Bitewing holder



Xios Endo holder



## Xios Intraoral Holder Systems Area-specific Rings & Arms

Easy, hygienic, and efficient solution for positioning digital intraoral imaging sensors

- To be used with all wired and wireless sensors
- Autoclavable Rings and Arms
- Single-use adhesive holder, eliminating any need for reprocessing
- Practical positioning for efficient workflow

Xios Anterior  
Ring & Arm



Xios Posterior  
Ring & Arm



Xios Bitewing  
Ring & Arm



## Xios Intraoral Holders for Seamless Intraoral Imaging

- Disposable holders – eliminating any need for reprocessing and save time
- Hygienic & Safe – A new sticky tab for every patient
- Consistent & Practical Positioning – Adhesive holders are designed for easy sensor placement, ensuring an efficient workflow
- Color-coded, area-specific holders simplify sensor positioning



Xios Sensor Sheaths


### Xios Sensor Holder Starter Kits

- 6411289 Xios XG Starter Kit S0-S1  
Contains: \*posterior ring (1), \*anterior ring (1), \*bitewing ring (1), \*posterior arm (1), \*bitewing arm (1), \*anterior arm (1), anterior holders blue (15), bitewing holders red (15), periapical yellow holders (15) endo green holders (50), Xios XG sheaths S0/S1
- 6411297 Xios XG Starter Kit S2  
Contains: \*posterior ring (1), \*anterior ring (1), \*bitewing ring (1), \*posterior arm (1), \*bitewing arm (1), \*anterior arm (1), anterior holders blue (15), bitewing holders red (15), periapical yellow holders (15) endo green holders (50), Xios XG sheaths S2 (50)

\* Products marked with (\*) are included as part of the Starter Kits and are not available for individual purchase. Please refer to the Starter Kit description for details on included components.

Xios Holders			Sensor Sheaths for Xios		
		Qty			Qty
6176528	Xios posterior sensor holder, yellow	100	6409960	Hygienic covers, Xios AE/XG, Size 0/1	300
6176510	Xios anterior sensor holder, blue	100	6409952	Hygienic covers, Xios AE/XG, Size 2	300
6176536	Xios bitewing holder, red	100			
6176544	Xios endo universal holder, green	50			
6176551	Xios endo sensor holder, grey	50			

## Xios Holders vs. Rinn® XCP Holders Color Guide



*Xios Anterior	*Xios Posterior	*Xios Bitewing	Xios Endo grey	*Xios Endo green
RINN Anterior	RINN Posterior	RINN Bitewing	RINN Endo	

\* Holders marked with (\*) are included as part of the Starter Kits. Please refer to the Starter Kit description for details on included components.

Universal & Fast Loading Design to fit all Digital Sensors.



## Rinn® Snap-A-Ray® DS Universal Sensor Holder

Fast and easy to use.  
For intraoral radiographs with bisecting technique.

- Sliding lock holds 1, 1.5, and 2 securely
- For anterior and posterior radiographs
- Autoclavable
- Snap-A-Ray® DS Endo features a cut out bite surface to accommodate endodontic files




Rinn® Snap-A-Ray® DS Universal Sensor Holder		Qty
550321	Snap-A-Ray® DS	1 pc.
550323	Snap-A-Ray® DS	3 pcs.
550421	Snap-A-Ray® DS Endo	1 pc.
550423	Snap-A-Ray® DS Endo	3 pcs.

Rinn® Snap-A-Ray® DS Positioning Arm and Ring Kit	
550322	Snap-A-Ray® DS Positioning Arm and Ring Kit Contains each: 2-in-1 aiming ring (1), Straight arm (1), Angled arm (1)

## Rinn® Eeze-Grip® Digital Sensor Holder

- Fits most sensor brands and sizes
- Great for all anterior and posterior radiographs, including the third-molar region
- Patented flexible grips designed for patient comfort
- Autoclavable



Rinn® Eeze-Grip® Digital Sensor Holder		Qty
550293	Eeze-Grip® Digital Sensor Holder	1 pc.
550294	Eeze-Grip® Digital Sensor Holder	3 pcs.

## Rinn® Uni-Grip® 360 Universal Sensor Holder


84% of sensor users found Uni-Grip® 360 to be faster than individual holders.<sup>1</sup>

- Designed to capture an entire full mouth series using only one holder
- Single-use holder accommodates any #1, #1.5 or #2 sized digital sensor
- Easy rotating tab for posteriors, anteriors and bitewings
- Perforated adhesive tabs for fast, easy removal
- Enhanced patient comfort due to compact size and round corners

Rinn® Uni-Grip® 360	
550052	Uni-Grip® 360 Kit Contains: 50 Uni-Grip® 360 Universal Sensor Holders, 50 Sensor Covers
550252	Uni-Grip® 360 Combo Pack Contains: 100 Uni-Grip® 360 Universal Sensor Holders, 100 Sensor Covers, 1 Uni-GripAR® Positioning Arm and Ring




Sliding the bitepiece along the channels enables you to take anterior, posterior and bitewing x-rays.



<sup>1</sup> Dentsply Sirona Internal User Study, 2012; N=24 clinicians.

# Less Shifting. More Security.


Designed to keep phosphor plates in place.



## Rinn® XCP-PSP Fit® Phosphor Plate Holder System

Less shifting and movement –  
minimizing the risk of errors and retakes

- Secure band envelopes work together with the biteblocks to keep the phosphor plate in place.
- Fits all major phosphor plate brands in sizes 0, 1, 2 and 3
- Protects phosphor plates
- Simplifies with color-coding assembly
- Enhances patient comfort with soft envelopes and rounded biteblocks
- Decreases clutter when used with Rinn® XCP-ORA® One Ring & Arm Positioning System



## Rinn® XCP-PSP Fit® Phosphor Plate Envelopes with Secure Band

Envelopes with Secure Band		Qty
552714	Phosphor Plate Envelopes, Size #2	100 pcs.
552710	Phosphor Plate Envelopes, Size #0	100 pcs.
552711	Phosphor Plate Envelopes, Size #1	100 pcs.
552712	Phosphor Plate Envelopes, Size #2	200 pcs.
552713	Phosphor Plate Envelopes, Size #3	100 pcs.



## Rinn® XCP-PSP Fit® Biteblocks



Anterior Thick		Fits Phosphor Plate Size	Qty	Bitewing Horizontal		Fits Phosphor Plate Size	Qty
552730	Anterior Thick	1	2 pcs.	552760	Bitewing Horizontal	0	4 pcs.
552731	Anterior Thick	2	2 pcs.	552761	Bitewing Horizontal	1	4 pcs.
<b>Anterior Thin</b>				552762	Bitewing Horizontal	2	4 pcs.
552740	Anterior Thin	0	4 pcs.	552763	Bitewing Horizontal	2	10 pcs.
552741	Anterior Thin	1	4 pcs.	552764	Bitewing Horizontal	3	4 pcs.
552742	Anterior Thin	2	4 pcs.	<b>Bitewing Vertical</b>			
552743	Anterior Thin	2	10 pcs.	552770	Bitewing Vertical	2	4 pcs.
<b>Posterior</b>				<b>Endodontic</b>			
552750	Posterior	0	4 pcs.	552781	Endo Horizontal	2	4 pcs.
552751	Posterior	1	4 pcs.	552790	Endo Vertical	1	4 pcs.
552752	Posterior	2	4 pcs.	552791	Endo Vertical	2	4 pcs.
552753	Posterior	2	10 pcs.				

## Rinn® XCP-ORA® One Ring & Arm Positioning System




## Rinn® XCP® Aiming Ring & Arm



Rinn® XCP-ORA® Ring & Arm		Qty	Rinn® XCP® Endo Arm & Ring		Qty
550771	XCP-ORA® Arm & Ring	1 pc.	550597	Endo Arm (No pins)	1 pc.
550772	XCP-ORA® Arm & Ring	2 pcs.	550598	Endo Ring (Green)	1 pc.
550773	XCP-ORA® Ring	1 pc.			
550774	XCP-ORA® Arm	1 pc.			

Angled for  
hard-to-reach areas



## Rinn® Snap-A-Ray® Xtra Film & Phosphor Plate Holder

### Film and Phosphor Plate Holder for Bisecting Technique

- Cushioned media grips hold film and phosphor plates secure!
- Cushioned bite areas designed to enhance patient comfort
- Media can be angled for hard-to-reach areas
- Autoclavable

Rinn® Snap-A-Ray® Xtra Film & Phosphor Plate Holder	Qty
540501 Snap-A-Ray® Xtra Film & Phosphor Plate Holder	1 pc.
540503 Snap-A-Ray® Xtra Film & Phosphor Plate Holder	3 pcs.

## Rinn® XCP® Holder System

Recommended for Film



Rinn® XCP® Kits	
542045	XCP® Endo Kit Contains: XCP® Positioning Arm Endo (1), XCP® Positioning Ring Endo (1), XCP® Biteblocks Endo (3)


XCP® Biteblocks		Qty
540861	Anterior	25 pcs.
540862	Posterior	25 pcs.
540867	Anterior	6 pcs.
540868	Posterior	6 pcs.
540863	Anterior and Posterior	10 pcs.
540929	Bitewings Size #2, Horizontal	6 pcs.
540931	Bitewings Size #3, Horizontal	6 pcs.
540935	Endo	3 pcs.

XCP® Positioning Arms		Qty
540857	Anterior	1 pc.
540858	Posterior	1 pc.
540927	Bitewings	1 pc.
550597	Endo (No pins)	1 pc.




XCP® Positioning Rings		Qty
540865	Anterior	1 pc.
540860	Posterior	1 pc.
540934	Bitewings	1 pc.
550598	Endo	1 pc.

## Rinn® XCP-ORA® Complete Kit for Film

Rinn® XCP-ORA®	
542071	XCP-ORA® Complete Kit for Film (no endo included) Contains: 72 Biteblocks – Anterior (20) Posterior (20) Bitewings Size #0 & #1, Horizontal (10) Bitewings Size #1, Vertical (6) Bitewings Size #2, Horizontal (10) Bitewings Size #2, Vertical (6) ORA® Ring & Arm (1) Positioning Poster (1)



# Rinn® User Guide

Media in Use	<p>Digital Sensors</p> 	<p>Phosphor Plates</p> 	<p>Film</p> 
General Use	 <p>Rinn® XCP-DS Fit®</p>	 <p>Rinn® XCP-PSP Fit™</p>	 <p>Rinn® XCP®</p>
Endodontic Procedures	 <p>Rinn® XCP-DS Fit® Endo Kit    Rinn® Snap-A-Ray® DS Endo</p>	 <p>Rinn® XCP-PSP Fit™ Endo Kit</p>	 <p>Rinn® XCP® Endo Kit</p>
Pediatric Procedures	 <p>Rinn® XCP-DS Fit® Size 1 Vertical Bitewing to be used with size 0 sensors for horizontal Bitewings</p>	 <p>Rinn® XCP-PSP Fit™ Size 0 &amp; Size 1 Kits</p>	
Single Use Option	 <p>Xios Holder Rinn® Uni-Grip® 360</p>	 <p>Rinn® Uni-Grip®</p>	
Bisecting Options (for challenging situations such as Narrow Jawbone and Tori)	 <p>Rinn® Snap-A-Ray® DS</p>	 <p>Rinn® Snap-A-Ray® Xtra</p>	

# Rinn® Color Guide

Anterior		
Posterior		
Bitewing		
Endo	